



Sterile alkaline cleaner for use in cleanrooms

## perform<sup>®</sup> sterile mucasol

### Our Plus

- ready-to-use alkaline cleaner for powerful cleaning
- fast, effective, alkaline (pH ca. 11.6)
- sterile (X-ray irradiated)
- ready to use in handy, sterile 500 ml and 1 L spray system
- colorant & perfume-free
- suitable for cleaning cytostatics and monoclonal antibodies

### Application areas

For cleaning surfaces and equipment/appliances in cleanrooms as well as surfaces made of glass, porcelain, plastic, rubber, stainless steel, brass and copper. Especially for areas of cleanliness classes A and B, in accordance with: Annex 1 to the EU Guide to Good Manufacturing Practice: Manufacture of Sterile Medicinal Products, March 2009. FDA Guidance for Industry: Sterile drug products produced by aseptic processing - current Good Manufacturing practice; U.S. Department of Health and Human Services, Food and Drug Administration, September 2004.

### Instructions for use

Apply undiluted to a cloth or objects and surfaces and wipe. After cleaning, rinse objects and surfaces thoroughly with water of at least drinking water quality to remove cleaning agent residues. Avoid contact with aluminum. After cleaning, subsequent disinfection with, for example, perform<sup>®</sup> sterile alcohol IPA\* is recommended.

Special notes:

The sterile spray system offers the following advantages:

- particularly handy and space-saving
- ergonomically shaped spray pump does not slip out of the gloved hand
- sterile even after opening

The bottles are sterile (with irradiation indicator) and individually double-wrapped.

\* Use disinfectants safely. Always read the label and product information before use.



# perform<sup>®</sup> sterile mucasol

## Product data

### Composition:

100 g solution contains the following active ingredients: < 5% phosphates (tripotassium orthophosphate)

## Special advice

Intended exclusively for use by skilled personnel. The product is unsuitable for cleaning alkali-sensitive materials (e.g. anodized aluminum and zinc). Do not mix with other products. Protect from frost, heat and sunlight. Store upright in the original container at room temperature. Dispose of completely empty packaging to the collection of recyclable materials. Keep container tightly closed.

## Information for order

Item	Delivery form	Item no.
perform <sup>®</sup> sterile mucasol 1 l bottle	6/Carton	on request
perform <sup>®</sup> sterile mucasol 500 ml bottle	10/Carton	on request

These products are not available in every country. For more information please contact our local subsidiary or distributor.

## Environmental information

schülke manufactures products economically and with advanced, safe and environmentally friendly production processes while at the same time maintaining our high quality standards.

## Expert opinion and information

Please visit our website for an overview of all available literature/reports on the product: [www.schuelke.com](http://www.schuelke.com).

For individual questions:

Application Department

Phone: +49 40 52100-666

E-Mail: [info@schuelke.com](mailto:info@schuelke.com)



Schülke & Mayr GmbH holds a Manufacturer's Authorisation according to sect 13 para 1 German Drug Law and Certificates of GMP Compliance for medicinal products.

schülke Headquarters  
Schülke & Mayr GmbH  
Robert-Koch-Str. 2  
22851 Norderstedt, Germany  
Phone +49 40 - 52100 - 0  
Fax +49 40 - 52100 - 318  
[www.schuelke.com](http://www.schuelke.com)  
[info@schuelke.com](mailto:info@schuelke.com)