

# PRODUCT INFORMATION

## Berner Micro-Spike

Spike with extra-long cannula for mini-vials



### Scope of application and properties

- + **Maximum protection:** The withdrawal cannula **Berner Micro-Spike** was developed for use with CMR drugs such as cytostatics, virostatics, anabolics etc. and offer maximum protection in combination with easy handling.
- + **Area of application:** The withdrawal cannula **Micro-Spike** is used for non-drip multiple extraction from supply bottles (vials) and for injections.
- + **Distinctive features of the Micro-Spike:** Especially long cannula (stainless steel needle tubing in accordance with ISO 9626 suitable for medical devices); no particle filter (liquid filter); 0.1 µm air PTFE vent filter; Latex-, PVC- and DEHP-free; ergonomic design for optimum handling. Protective cap prevents contamination.
- + **Before use:** Inspect sterile blister packaging for any damage! Do not use if blister pack is damaged. We recommend disinfecting the blister pack using a sterilizing wipe before use. Single use product.
- + **Application instructions:** When injecting liquid do not hold the vial upside down, as it clogs the PTFE membrane and prevents pressure equalization.
- + **Disposal:** For use in connection with CMR<sup>1)</sup> drugs: Waste requiring supervision (European waste code: 18 01 04 in as per 2000/532/EC); in case of heavy contamination, waste requiring special supervision (European waste code: 18 01 08\* as per 2000/532/EC); separate collection and disposal!

<sup>1)</sup>: Carcinogenic mutagenic reproductive toxic

### Version

<b>Order No.</b>	<b>126215</b>	<b>PU</b>	<b>50 pieces</b>
------------------	---------------	-----------	------------------

- + No particle filter (liquid filter)
- + 0.1 m bacteria-free PTFE aeration and deaeration filter (aerosol filtration efficiency: 99.999997%)
- + long stainless steel cannula (according to EN ISO 9626:1995), Ø 1.5 x 45 mm
- + Latex-, silicone oil-, PVC- and PHT/DEHP-free,
- + Housing: transparent MABS
- + Low dead volume: 0.115 ml
- + Integrated protective cap as contamination protection

### CE marking & notified body

CE-marking in accordance with medical devices directive 93/42/EEC. Notified body "0482"; MEDCERT Certification and Notified Test Body for the Medicine GmbH, Pilatuspool 2, 20355 Hamburg

### Quality management system

Our **quality management system** is **tested and certified** by the TÜV Product Service GmbH in accordance with DIN EN ISO 9001:2015. Regular **audits** and **production site inspections** guarantee the quality of our products.

### Storage and transport conditions

Dark (protect against UV- and sun light), cool (+15 to +25°C), dry (relative humidity 30% - 60%)

### Shelf life

5 years from the date of sterilisation

---

## Manufacturer

ROWEMED AG; Juri-Gagarin-Ring 4; D-19370 Parchim, Germany

---

## Distributor

Berner International GmbH, Werner-von-Siemens-Str. 19, D-25337 Elmshorn, Germany

---

## Instruction manual

