# SpillKit XP

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## 1. User information

Read the user information carefully and keep it in the SpillKit. In accordance with the applicable legal guidelines, e.g. for example the German Employers' Liability Insurance Association for Health Service and Welfare Work (BGW) and the „Quality standard for pharmaceutical oncological service (Quapos) all personnel who are involved in handling cytostatics should be taught how to use the emergency cleaning kit at least once a year. Ignoring the user information can impair the performance of the SpillKit components and be harmful to the health of the user and the environment.

## 2. General information

### 2.1 Manufacturer Berner International GmbH

### 2.2 Quality Assurance

Tested and certified in accordance with EN ISO 9001:2015 by the TUV Product Service GmbH (certificate see [www.Berner-safety.de](http://www.Berner-safety.de)).

### 2.3 Copyright and Industrial Property Rights

The product group symbol is a registered and protected trademark. Infringement of industrial property rights will be prosecuted. Berner International GmbH assumes no warranty or liability for Berner SpillKit components that are employed beyond the scope of their intended use. For translations into foreign languages, the German version of the user information is binding. The content can be changed at any time and without notice. This user information is the property of Berner International GmbH and may not be forwarded to third parties, copied, photocopied or used for promotional purposes, either wholly or in part, without written consent. Subject to technical and production-related changes. © Copyright Berner International GmbH

## 3. Description

The emergency cleaning kit is suitable for the safe disposal of small quantities of spilt CMR drugs (e.g. cytostatic drugs). The Berner SpillKit contains a complete set of personal protective equipment, a warning sign and a pen for marking out the accident site, cleaning utensils for dry and liquid contamination as well as broken glass, waste disposal utensils and an accident report form (only relevant for Germany). The Berner SpillKit contains the components listed in section 4.2 of QuapoS and therefore complies with German Society of Oncology Pharmacy (DGOP) and the European Society of Oncology Pharmacy (ESOP) recommendations for a SpillKit.

### 3.1 Contents Berner SpillKit

Order No.	Item Description	Quantity
4104	Berner SpillKit XP, Size M	1
4112	Berner SpillKit XP, Size L	1
4115	Berner SpillKit XP, Size XXL	1
4124	Berner SpillKit XP latexfree, Size M	1
4133	Berner SpillKit XP latexfree, Size L	1
4135	Berner SpillKit XP latexfree, Size XXL	1

Item Description	Quantity / Pair		Quantity / Pair
<b>Personal Protective Equipment (PPE)</b>		<b>Tools</b>	
Cytostatic protective overall*	1	Scoop	1
Cytostatic protective gloves**	1	Sweeper	1
Protective gloves**	1	Tongs	1
Protective overboots	1	Bottle, distilled water (500 ml)	1
FFP3 respirator mask	1	Waste bag, transparent	1
Protective gown	1	Special waste bag, white	1
<b>Tools</b>		Cable ties	2
Marker pen	1	LDPE Minigrip bag	1
Warning sign	1	User guide	1
ChemoSorb Pad	1	Accident report (only for Germany)	1
ISYSOFT Wipes	6	Cardboard box	1
* Size M - 4104 / 4124; Size L- 4112 / 4135; Size XXL- 4115 / 4135 **Latex or Latexfree			

### 3.2 Usage

The emergency and cleaning kit must only be used once to dispose of spilled cytostatic drugs. Do not reuse any parts of the SpillKit contents!

### 3.3 Disposal

After use, all Berner SpillKit components should be treated as waste requiring special supervision (waste code: 18 01 08\* in accordance with 2000/532/EC). Use the waste disposal utensils provided. Only approved specialist waste disposal companies may collect, transport and dispose of used kits.

### 3.4 Storage and Transport

Dark (away from direct UV light and sunlight); cool (+5 to +40°C); dry; no contact with pointed and/or sharp objects. Away from equipment or installations that can produce ozone (e.g. through mercury vapour lamps, high-voltage equipment, etc.); avoid direct contact with metals, e.g. copper, manganese, magnesium and iron; avoid contact with oil-based antiseptic phenols and their derivatives, fats, petrolatum, petroleum, paraffin or other similar compounds. Breakable.

### 3.5 Durability

The components of the SpillKit are subject to aging. This process known as degradation can lead to the loss of protective functions or general effectiveness. SpillKits which have exceeded their use-by date should therefore not be used for the clean up of „real“ accidents/spills. We recommend using these SpillKits for training and instruction purposes of personnel.

### Shelf life

**Complete Berner SpillKit: 3 years from date of manufacture.**

### 3.6 LOT-Number

LOT-No.: 4 Digits, allocated consecutively.

## 4. Instructions for use

The EC Directive 2000/54/EC and the relevant national regulations (e.g. Germany: ChemG, GefStoffV, TRGS 525) regulate the necessary behaviour in situations where employees are exposed to unusually high concentrations of carcinogenic hazardous substances. Appropriate precautions must be taken in this regard. The employer has the duty and responsibility to carry out a hazard analysis in advance, to draw up the operating instructions and to instruct all relevant employees practically at least once a year. This must be documented in writing. Have a Berner SpillKit ready and attach a short user manual. Analyse possible hazardous situations in advance, e.g. in the context of transport, delivery, storage, manufacture, application, cleaning, disposal, etc.

**In the event of an accident resulting in cytostatic contamination, please follow the instructions below:**

#### Step 1: Keep calm!

Always deal with personal contamination before material contamination. In the event of contact with skin, rinse immediately in plenty of cold running water. Do not use soap. In the event of contact with the eyes, rinse immediately with plenty of water or isotonic saline solution for at least 10 minutes. Remove contaminated pieces of clothing immediately. If possible, ask someone else for help. Please note that decontamination should only be performed by trained staff.

#### Step 2: Follow the brief directions for use / user guide!

Follow the brief directions for use/user guide, i.e. the quick user guide should be kept near the workplace and clearly visible. Carry out the individual steps, in accordance with the yearly training session.

#### Step 3: Don PPE!

First protect those involved. Put on the appropriate PPE and ensure that the person assisting you also has sufficient protection.

**Protective clothing must be put on in the correct order:**

- 1. Put on the respiratory mask.** Hold the mask in your hand with the straps hanging down(1). Guide the first strap over your head as far as the nape of your neck . Place the mask under your chin. Guide the the second strap over the back of your head (2). Adjust the nosepiece to fit your nose (3). Cover the half-mask with both hands and breathe out strongly (4). The mask material will move up and down. If air escapes at the edge of the mask, move the mask into the correct position and adjust the nosepiece to fit your nose. Please note that facial hair and scars can impair the fit of the mask.
- 2. Don the safety goggles.** Use the head strap to adjust them and ensure a secure fit. Safety goggles can be worn over corrective glasses.
- 3. Put on the protective overall** for use with cytostatic drugs: Undo the zip and put your arms and legs inside the overall.
- 4. Put on the blue or buff-coloured protective gloves for use with cytostatic drugs** and pull the cuffs over the ends of the overall sleeves.
- 5. Then, put on the protective overboots.**
- 6. Lastly, don the green-coloured protective gloves** over the blue or buff coloured gloves.

#### **Step 4: Mark clearly and close off the accident site!**

Secure the site of the accident to prevent the spread of contamination (this applies particularly to freely accessible areas, e.g. stations, corridors, etc.). Liberally mark out the accident site using the pen provided and close off the affected area with the warning sign. Take steps to prevent turbulence: switch off the ventilation system, close windows and doors and stop people using the area. In the event of a spillage inside a safety cabinet, switch on the ventilation and leave it running.

#### **Step 5: Dispose of the spillage and clean as directed!**

Please note that the gloves do not usually provide adequate protection against cuts. The tongs or sweeper should therefore be used for glass and sharp fragments. Liquid contamination is absorbed using ChemoSorb pads and cloths: spread the pad out over the spill liquid and it will soon turn into a jelly-like mass. For contamination in powder form, use the distilled water and ISYSOFT cloths: place the dampened cloths over the spillage. Always clean from the outside in. Do not use a hand brush or other similar tool as this produces turbulence

#### **Step 6: Place the contamination in the transparent waste bag!**

Pick up glass and sharp fragments with the tongs and place in the box. Use the spatula and scoop for pads and cloths. Always carry out several stages of cleaning (basic cleaning-> follow-up cleaning, see section 6). First place the contamination in the box. Then put the box, tools and cleaning utensils in the transparent waste bag and seal with a cable tie.

#### **Step 7: PPE and transparent waste bag in white special waste disposal bag!**

Place the transparent waste bag in the white special waste bag. Take off the protective clothing and place it in the white special waste bag. This must be done in the correct order: 1. Overboots 2. Overgloves 3. Overall 4. Cytostatic gloves 5. Safety goggles 6. Respiratory mask. Now close the white special waste bag with the second cable tie. Follow the disposal instructions in section 3.3.

#### **Step 8: Shower!**

Shower thoroughly.

#### **Step 9: GP**

If necessary see a doctor, especially in case of personal contamination.

#### **Step 10: Accident Report**

Write a precise accident report. Inform the employee responsible for health and safety at work and the company doctor about the accident. After an accident, hold a meeting and/or further training if necessary.

**In all cases, ensure you obtain a new Berner SpillKit.**

## **5. Product properties**

### **5.1. activebreath 2 coverall and protective overboots**

#### **Properties coverall**

Maximum protection and wearing comfort: Type-tested and certificated as complex PPE category III; Chemical protective equipment type 5/6; optimum personal and product protection; comfortable to wear (breathable); material is low in lint and particles, anti-statically treated. High neckline, zipper with overlapping flap, thumb loop, elasticated sleeves and ankles, hood. Application range: Protective overall suitable for the handling of CMR pharmaceuticals (e.g. cytostatics, virostatics). Protective barrier: Reliable barrier against numerous water-based inorganic chemicals at low concentration and particles larger than 1 µm. (category III). Protection against infectious agents. Protection against antistatic particles (EN 1149-5). Protective properties: Protection against all CMR drugs or chemicals cannot be guaranteed! In case of exposure to biohazards, which do not correspond to the degree of imperviousness of the protective clothing, bio-contamination of the wearer is possible. Keep away from naked flames and heat sources.

#### **Properties overboots**

Maximum protection and comfort: Type-tested and certified as complex PPE of the highest category III; chemical protective clothing type PB [6]-B, partial body protection. Optimal personal and product protection; liquid-impermeable in the coated foot area, pleasant wearing comfort; overboots with non-slip sole, leg closure with elastic band; additional banding for fixation; sterile and non-sterile design Application range: The protective overboots protect against contact with CMR drugs. They provide protection against particles, liquids and aerosols type PB [6]-B, with limited protection against liquid chemicals Protection barrier: Liquid impermeable coating with a high barrier function of the coated material against bacteria and viruses. Protective properties: Protection against penetration of chemicals and biologically contaminated liquids, limited spray tightness according to EN 13034, tested permeation protection according to the permeation list given here. No guarantee for CMR drugs or chemicals not listed. In case of exposure to biological agents and chemicals, which do not correspond to the level of protection of the protective clothing, permeation of the protective overboots type PB [6]-B may occur.

### Types protective coverall

SpillKit XP	Size	Chest (cm)	Height (cm)	Colour
4104 / 4124	M	92 – 100	170 – 176	white
4112 / 4133	L	100 – 108	176 – 182	white
4115 / 4135	XXL	116 – 124	188 – 194	white

### Types protective overboots

Size	unisize	Colour	white
------	---------	--------	-------

### Material properties activebreath 2 coverall and overboots

Material	Non-woven polypropylene
Material weight	69 g/m <sup>2</sup>
Coverall seams	untaped

### Protection against mechanical hazards

Requirement	Performance level
Abrasion resistance (1-6) to EN 530:2010	2
Puncture resistance (1-5) to EN 863:1995	1
Seam resistance (1-5) to ISO 13935-2:2014	3
Trapezoidal tear strength (1-6) to EN ISO 9073-4	6

### Protection against chemical hazards

Permeation <sup>1)</sup> tested to EN 14325. Determined for the following substances:	
5-Fluorouracil, 10 mg/ml	Class 4
Cyclophosphamide, 20 mg/ml	Class 4
Isopropanol, 70%	Class 2
NaOH, 10%	Class 4
Paclitaxel, 6 mg/ml	Class 4
Sulphuric acid, 30%	Class 3
Thiotepa, 10 mg/ml	Class 4
Vincristine, 20mg/ml	Class 4
<sup>1)</sup> : Movement of a chemical through material on a molecular level. The performance class does not reflect the actual duration of protection at the workplace!	

### Care instructions

Do not wash, do not iron, do not tumble dry, do not dry clean, single use only.

### CE-Marking

CE marking acc. of the EU PPE regulation 2016/425 for complex PPE category III. Protective overall: Type test carried out according to EN 13982-1: 2004, EN 13034: 2005 + A1: 2009 and EN 1149-5: 2008; Notified body 0624. Protective overboots: Tested in accordance with EN 13034: 2005 + A1: 2009, EN 14126: 2003; Notified body 0120.

## 5.2. Protective gloves Berner Manu L and Manu Prene XP

### Properties

Maximum protection and comfort: type tested and certified as complex highest level PPE of the category III; good grip, good tactile sensitivity. Area of application: Protective gloves for handling CMR drugs (e.g. cytostatics, virostatics) and biological agents (e.g. bacteria and viruses). Protective capacity: No guarantee for all CMR drugs or chemicals! In case of exposure to biological hazards, which do not correspond to the degree of imperviousness of the protective clothing, bio-contamination of the wearer is possible. Before use: Check for any damage! Do not use damaged gloves!

### Types

SpillKit XP	4104	4112	4115	4124	4133	4135
Size	7.5 / M	8.5 / L	9 / XL	7.5 / M	8.5 / L	9 / XL
Protective gloves	Manu L			Manu Prene XP		

### Material properties

#### 5.2.1 Protective glove Berner Manu L

Natural latex (latex allergen < 0.5 µg/g); low protein (protein = 17 µg/g); powder free in accordance with TRGS 540; latex and carbamates can trigger allergies. AQL: 0.65

### Material thickness

Measurement points	Material thickness d (measured twice)
Finger, 15 mm from the end of the tip	≥ 0,96 mm
Middle of the palm	≥ 0,86 mm
Shaft, 25 mm from the end of the shaft	≥ 0,48 mm

**Protection against chemical hazards; permeation according to EN ISO 374-1:2016; test method EN 16523-1:2015 (\*tests according to EN 374-3); degradation 374-4:2013. Breakthrough times [min] / performance classes(1-6) were determined for the following chemicals:**

Chemical	Breakthrough time [min]	Performance class (1-6)	Degradation
40% Sodium hydroxide (K)	> 480	6	-6.5%
37% Formaldehyde (T)	> 480	6	-3.1%
30% Hydrogen peroxide (P)	> 480	6	-12.5%
Bleomycin 3 mg/ml*	> 180	4	-
Carboplatin 10 mg/ml*	> 180	4	-
Carmustine 4 mg/ml*	> 60	3	-
Isopropanol 70%*	> 30	2	-
Isopropanol 70% + Carmustine 4 mg/ml*	> 120	4	-
Cisplatin 50 mg/ml*	> 120	4	-
Cyclophosphamide monohydrate 20 mg/ml*	> 480	6	-
Doxorubicin hydrochloride 1 mg/ml*	> 120	4	-
Daunorubicin hydrochloride 1,5 mg/ml*	> 60	3	-
5-Fluorouracil 1,5 mg/ml*	> 60	3	-
Methotrexate 2 mg/ml*	> 120	4	-
Mitomycin 1 mg/ml*	> 180	4	-
Vinblastine 1 mg/ml*	> 180	4	-
Vincristine 1 mg/ml*	> 120	4	-
Sulphuric acid 40%*	> 480	6	-
Sulphuric acid 96%*	> 30	2	-
Sodium hydroxide 10%*	> 480	6	-
Sodium hydroxide 30%*	> 480	6	-

Glutaraldehyde 5%*	> 480	6	-
Diethylamine (undiluted)*	> 30	2	-

#### Protection against biological hazards; penetration in accordance to EN 374-2:2014:

Tested according to EN 374-2:2014- test conditions fulfilled. Air leak test only

#### Resistance against viruses, bacteria & fungi; tested against blood-borne pathogens according to EN 374-5:2016

Test conditions fulfilled. In accordance with current knowledge, it should be assumed that meeting the penetration requirements provides effective protection from microbiological hazards.

#### CE-Marking

CE marking according to EU PPE Regulation 2016/425 for complex category III PPE. Type examination performed based on EN ISO 374-1:2016 Type B; EN 16523-1:2015, EN 374-2:2014, EN 374-4:2013; EN 374-5:2016; EN 420:2003+A1:2009-11; EC type examination certificate No. CE 710563;

**Notified Body 2797:** BSI Group, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes, MK5 8PP, UK

#### 5.2.2 Protective glove Berner Manu N

Latex-free protective glove made from Nitrile; powder-free in accordance with TRGS 540

##### Material thickness

Measurement points	Material thickness d (measured twice)
Finger, 15 mm from the end of the tip	≥ 0,28 mm
Middle of the palm	≥ 0,18 mm

**Protection against chemical hazards; permeation according to EN 374-1:2016; test method EN 16523-1:2015 (\*tests according to EN 374-3); Degradation 374-4:2013. Breakthrough times [min] / performance classes(1-6) were determined for the following chemicals:**

Chemical	Breakthrough [min]	Performance class (1-6)	Degradation
40% Sodiumhydroxide (K)	> 480	6	-34.0%
37% Formaldehyde (T)	> 241	5	8.0%
30% Hydrogenperoxide (P)	> 121	4	-0.8%
5-Fluorouracil 10 mg/ml*	>480	6	
Carmustine, 3.300 ppm*	90	3	
Cisplatin, 1.000 ppm*	>480	6	
Cyclophosphamide monohydrate, 20.000 ppm*	>480	6	
Doxorubicin hydrochloride, 2.000 ppm*	>480	6	
Etoposide, 20.000 ppm*	>480	6	

#### Protection against biological hazards; penetration in accordance to EN 374-2:2014

Tested according to EN 374-2:2014- test conditions fulfilled.

#### Resistance against viruses, bacteria & fungi; tested against blood-borne pathogens according to EN 374-5:2016

Test conditions fulfilled. In accordance with current knowledge, it should be assumed that meeting the penetration requirements provides effective protection from microbiological hazards.

#### CE-Marking

CE marking according to the EU PPE regulation (EU) 2016/425 for complex PPE of category III. Type examination performed based on EN ISO 374-1:2016 Type B; EN 16523-1:2015, EN 374-2:2014, EN 374-4:2013; EN 374-5:2016; EN 420:2003+A1:2009; EC Type Examination Certificate No. 2777/10922-02/E02-01.

Notified body: SATRA Technology Europe Ltd, Bracetown Business Park, Clonee, D15YN29P, Republic of Ireland.

#### 5.2.3 Protective glove uvex profastrong NF33

##### Properties

Chemical protective glove type A, anatomically shaped, high wearing comfort, good grip and sensitivity; AQL= 1.5; very good mechanical capacity. Protective capacity: No guarantee against all chemicals! Area of application: To be worn over the protective gloves for use with cytostatic drugs. Before use: Check for any damage! Do not use damaged gloves!

## Types

SpillKit XP	4104	4112	4115	4124	4133	4135
Size	9	9	10	9	9	10
Protective glove	Nitrile, latexfree					
Colour	green					

## Material properties

Glove length	330 mm
Material thickness	0.38 mm
Material	Nitrile, inside flock-lining

## Protection from mechanical hazards, in accordance to EN 388:2016. Coding as follows:

Requirement	Performance class (0-4)
Abrasion (0-4)	4
Cut (0-5)	1
Tear (0-4)	0
Puncture (0-4)	1

## Protection from chemical hazards; permeation in accordance with EN 374-1:2016, Type A:

Chemicals	Breakthrough time [min]	Performance class(1-6)	Degradation (%)
Methanol (A)	>30	2	48.4
n-heptane(I)	>480	6	18.7
Sodium hydroxide 40% (K)	>480	6	-9.1
Sulphuric acid 96% (L)	>120	4	72.7
Ammonium hydroxide 25% (O)	>480	5	3.1
Formaldehyde 37% (T)	>480	5	35.4

## Resistance to viruses, bacteria & fungi: Tested according to EN ISO 374-5:2016.

Tests passed.

## CE-Marking

According to the PPE Regulation (EU) 2016/425 for complex category III PPE based on EN 420:2003+A1:2009, EN 374-5:2016, EN 374-1:2016, EN 388:2016. Type examination and certification by the notified body no. 2777.

## Notified body

2777: Satra Technology Europe Ltd, Bracetown Business Park, Clonee, Dublin 15, Dublin, Ireland; quality assurance and monitoring by SGS Fimko Oy, Notified body 0598, Takomotie 8, FI-00380 Helsinki.

## 5.3 FFP3 respiratory mask

### Properties

Half-mask with high filtering capacity, tight-fitting and leakproof. Its soft material, special shape and nose clip ensure it can be easily adjusted to fit the face. Minimal respiratory resistance and the exhalation valve make it extremely comfortable to wear. The threaded elastic strap makes it easy to put on and take off. Area of application: protection from solid and liquid particles, e.g. for handling CMR drugs. Area of application: Protection from solid and liquid particles, e.g. for handling CMR drugs. Before use: Check for damage. Do not use damaged masks. Limitation: laboratory test results can differ from values obtained in practice considerably and, as a result, the useful life of the mask may increase or decrease. Facial hair and scars can impair the fit of the mask and therefore the level of protection.

## Types

Size	Universal
------	-----------

## Material properties

Latex-free	complete mask
------------	---------------

## Degree of particle filtration

Degree of particle filtration	Minimum degree of filtration
Sodium chloride	> 99%
Paraffin oil	> 99%
Dolomite dust test for FFP3- respiratory masks (Dolomite DRB4/15b testaerosol)	passed
Gases	not suitable

## CE-Marking

In accordance with PPE Regulation (EU) 2016/425 for complex PPE category III, based on EN 149. Type-testing, certification and control measures by the notified body 0426.

## Notified body 0426

Italcert S.r.L.; v.Le Sarca, 336, 20126 Milano, Italy

## 5.4 Safety goggles

### Properties

Full vision goggles; non-fogging and can be worn over corrective glasses. The elastic strap ensures a firm and fit. Area of application: Protection from mechanical risks, liquids and coarse dust. Protection capacity: Does not provide protection from high-speed particles, laser beams, temperatures of over 55°C, arcing faults, molten metals. Before use: Check for any damage! Do not use goggles if label is absent or illegible.

### Version

Size	Universal
Colour	Transparent

## Marking of the frame and lens

W 166 DIN CE 0196	Protection against mechanical risks, liquids, large particles
W 1N DIN CE 0196	Non-fogging, colourless plastic lens

## CE-Marking

In accordance with the PPE regulation (EU) 2016/425 for complex PPE category II based on EN 166. Type-tested and certified by the notified body 0196. Manufacturer UVEX Arbeitsschutz GmbH, Fürth, Germany.

## Notified body 0196

DIN CERTCO, Gartenstraße 133, D-73430 Aalen, Germany.

## 5.5 Further components of the Berner SpillKit XP

ChemoSorb-Pads	For absorption and binding of liquids
ISYSOFT Wipes	Extremely absorbent cloths, especially suitable for absorbing spilt liquids.
Marker pen	Colour black
Warning sign	Cardboard sign, Size: 210 x 150 mm (upright), Warning instructions in German & English.
Scoop	Material: Plastic, size: 135 x 100 mm, colour: dark blue, nitro-resistant, beaded back, untoothed.
Sweeper	Material: Plastic, size 135 x 100 mm
Tongs	Wood, no chemicals in manufacturing process.
Bottle, distilled water	Plastic bottle, 500ml distilled water.
Waste bag	Material: Polyethylene, size 576 x 1000 mm, thickness: 50 µm, colour: transparent, weight: 46 g/m <sup>2</sup> ± 5g.
Special waste bag	Polyethylene with print: „Zytotoxische und zytostatische Arzneimittel / Cytotoxic and cytostatic drugs / Farmaci Citotossici e Citostatici / Medicamentos Citotóxicos y citostáticos / Médicaments dytotoxiques et cytostatiques“, colour white, weight: 127 g ± 5g
Cable tie	Plastic, size 200 x 3,6 mm
LDPE-Minigrip-bag	Material, polyethylene, size: 300 x 400 x 0,05 mm, colour: transparent, resealable.

## 6. Cleaning

Cytostatic contamination must be cleaned thoroughly several times. In the main, testing focused on cleaning contaminated primary packaging although the results also provide a starting point for the problems caused by accident-related cytostatic contamination. The following table contains the results of 10 individual swab tests. Different cloths and solvents were tested on a glass plate:

Solvent used for swab test		Proportion of substance removed (Mean recovery rate in eluate)		
Stage 1	Stage 2	Fluorouracil	Methotrexate	Etoposide
Ethanol	Isopropanol	85%	59%	94%
0,05-0,1 M Caustic soda	Isopropanol	99%	95%	92%

Preference should therefore be given to two-stage cleaning (using a combination of caustic soda and isopropanol) over one-stage cleaning. The more frequent the cleaning procedure (at least 3-4 times), the more effective it is (95-99%).

**Removal of pen markings:** The markings usually disappear completely during the multi-stage cleaning process. Use isopropyl alcohol, for example, for removal.

# Use of SpillKits for the disposal of spilt cytotoxins

1

Keep calm!

2

Follow instructions!

3

Wear suitable PPE Cat. III !

4


Seal off and mark scene of accident!

5


Observe type of cleanup!

6

2-step cleaning!

7

Waste in blue disposal bag!

8

Waste bag and PPE in hazardous waste bag!

9

Shower thoroughly!

10

See a doctor!

11

Document incident!

Berner International GmbH  
Werner-von-Siemens-Str. 19  
25337 Elmshorn  
Germany

Tel.: +49 4121 4356 - 0  
Email [info@berner-safety.de](mailto:info@berner-safety.de)

[www.berner-safety.de](http://www.berner-safety.de)