

PRODUCT INFORMATION

Mop covers for cleaning tool II

Application and properties

Mop covers for the separately available isolator and cleaning tool for efficient and uncomplicated cleaning of hard-to-reach areas. For easier wall, floor and ceiling cleaning inside a safety cabinet. Available sterile packed, or pre-soaked with isopropanol.

Material properties cleaning cover

Material	microfibre polyester, multi-layer, bound pockets
Dimensions	ca. 20 x 7 cm
Version	beta sterilised in double bag, 2 pieces per pack



To order

	Order no.
Mop covers, sterile (2 pieces per pack):	
Double packed for cleanroom use	121742
Pre-saturated with mit Isopropanol	100100
Suitable for:	
Cleaning tool II	200662