

PRODUCT INFORMATION

Cleaning tool for isolators & safety cabinets

Applications and Features

- + Tool for isolators / glove boxes and safety cabinets for cleaning difficult to reach places efficiently and effectively. Walls, floors and ceilings are quick and easy to keep clean.
- + Made from stainless steel, electro-polished, with a double plastic hinge.
- + Short stainless steel for ease of handling and use.
- + Autoclavable at 121° C.
- + Resistant to disinfectants

Material Properties Cleaning Tool II

Material	Stainless steel, electro-polished, Hinge made from autoclavable plastic
Length of handle	ca. 580 mm
Total length	ca. 650 mm
Dimensions	ca. 190 x 50 mm

Material Properties Mop Cover

Material	Micro fibres polyester, several layers, fitted
Dimensions	ca. 20 x 7 cm
Beta sterilised, double bagged, 2 wipes per pack	



To order

	Order No.
Cleaning tool II	200662
Mop covers, sterile (2 per pack):	
Dry cloths	121742
Pre-saturated with Isopropyl alcohol	100100