

| | 126217 | 126218 | 126215 | 126213 | 126212 | 126211 | 126216 |
|------------------------------|------------------------|------------------|---------------------------------|-----------------|---------------|------------|--|
| Bezeichnung | Ultraspike II | Ultraspike II oF | MicroSpike | UltraSpike SKoF | UltraSpike SK | UltraSpike | Spike II |
| Anschluß | Luer-Lock | Luer-Lock | Luer-Lock | Luer-Lock | Luer-Lock | Luer-Lock | Luer-Lock |
| Metalldorn | | | ✓ | ✓ | ✓ | | |
| Kunststoffdorn | ✓ | ✓ | | | | ✓ | ✓ |
| Halkey-Roberts-Ventil | ✓ | ✓ | | ✓ | ✓ | ✓ | |
| Verschlusskappe Snap Fit | | | ✓ | | | | ✓ |
| Klammersystem | ✓ | | | | | | ✓ |
| Flüssigkeits-/Partikelfilter | | | | | ✓ | ✓ | |
| Be- und Entlüftungsfilter | ✓ | | ✓ | ✓ | ✓ | ✓ | ✓ |
| Füllvolumen | 0,2 ml | 0,25 ml | 0,2 ml | 0,25 ml | 0,25 ml | 0,25 ml | 0,13 ml |
| Restvolumen | 0,053 ml | 0,053 ml | 0,115 ml | 0,059 ml | 0,044 ml | 0,075 ml | 0,034 ml |
| Eignung | Viskose Substanzen | ✓ | ✓ | ✓ | ✓ | | ✓ |
| | Monoklonale Antikörper | ✓ | ✓ | ✓ | ✓ | | ✓ |
| | | | Metalldorn 1,5 mm DM X 45 mm | | | | Speziell für 20 mm Verschlusskappen |